

# ElectroTherapy Massager

## MV-950

### Relaxation, Anywhere, Anytime

Warranty Questions: [support@gurinproducts.com](mailto:support@gurinproducts.com)


#### Included in this standard package:

- Device x 1pc
- Large pad x 2pcs
- Medium pad x 4pcs
- Small pad x 2pcs
- Electric wire x 2pcs
- USB cable & AC adapter x 1set
- Instruction manual x 1pc
- Collection board for pad and wire x 1pc
- Bag x 1pc

#### Features

- 2 output with 16 modes
- 20 level intensity
- 20 minutes default timer, can adjust from 10-60 minutes
- Rechargeable lithium battery
- Battery indicator to mark when battery needs to recharge
- AC adapter 110-240v, USB Charger 5v

#### Illustrations of the device and accessories: (Model: MV-950)



- |                         |                          |
|-------------------------|--------------------------|
| 1. Two output connector | 8. Time adjustment       |
| 2. USB port             | 9. Pause switch          |
| 3. LCD screen           | 10. Pads- sticking board |
| 4. Battery indicator    | 11. Charger              |
| 5. On/off switch        | 12. USB cable            |
| 6. Intensity adjustment | 13. Massage pads         |
| 7. Mode selector        | 14. Pad wires            |

#### Modes

There are sixteen modes to choose from, there are:

- 1) Kneading/Knocking (4 Modes)
- 2) Acupuncture (3 Modes)
- 3) Cupping (3 Modes)
- 4) Tapping (2 Modes)
- 5) Shiatsu (2 Modes)
- 6) Low Intensity (2 Modes)

*Press the mode button to select mode. Once mode is selected, press the center button to choose the waves within that selection. Then slowly increase intensity to desire by pressing + and - buttons.*

#### Operating instructions

1. Connect the pads to the wire by snapping them on, then remove the protective film.

*Note: adding a few drops of water on the pads will increase their longevity.*

2. **Place two pads on the same area where the massage sensations are desired.** Plug pad wire into one of the ports on the device.

For additional relief use the second port and repeat procedure to put on another two pads.

Make sure **both** pads are on the skin and are not overlapping. Make sure your skin is free from any dirt, oil or lotion.

3. Turn the device on by sliding the on/off switch on the top of the unit. The unit will be now activated.

4. Press the mode button to try different mode. Then press the center button to choose the waves within that selection.

Remember the strength will default to the lowest setting by each mode selection.

*Note: Center button can also be used as a lock key.*

Press and hold center button for 3 second to lock the screen, and 3 second again to unlock the screen.

5. Once you select a mode, gradually increase the intensity by pressing the + button; decrease it by pressing - button.

6. Recommended application is 10-60 minutes. Increases by 10 minutes by each setting. The time remaining is automatically counted down. Then the unit will automatically turn off after time finishes.

## Precautions

- Do not apply the massage pads near the heart, on the head, above the neck, in the pubic region, over scarred areas, on the throat or over the mouth.
- Avoid touching the pads when the unit is on.
- Pregnant women must avoid using this massager.
- Do not use this massager if you have a pacemaker or other life support equipment such as an artificial lung or respirator.
- Do not use while driving, operating machinery, or during any activity in which involuntary muscle contractions may put the user at undue risk or injury.
- Before pads are applied to a new area, make sure device is turned off.

- Do not use if you have the following medical conditions:

Acute disease, trauma, or a surgical procedure in past six months

Cardiac dysfunction

Epilepsy

High blood pressure

High fever

Cancer

Infectious disease

Malignant tumors

Recent surgery

Skin problems

Pregnancy

Undiagnosed pain syndromes

Abdominal or inguinal hernia

Painful and/or atrophied muscles

Limited range of motion in skeletal joints

- Keep out of the reach of children.
- The unit is not intended for the application of any medical condition or disease nor is it intended for physiotherapy or muscle rehabilitation. It is contraindicated for use on any muscle that is injured or diseased.

## Maintenance and care

- When you are finished using the device, turn it off and take off the pads. Store the unit in a cool, dry place until next use.
- After usage, please cover the pads with the protective film and put them back in their zip lock bags. Applying a few drops water before and after usage will increase the longevity of the pads.
- Charge the unit at least every 3 months in the event it the device is not used regularly. You can also use the unit with the wall charger or directly from your computer using the USB cable.

## Troubleshooting

Problem	Cause	Solution
One pad feels stronger than the other.	This is normal. Different areas of your body will react differently.	Nothing needs to be done. Make sure the pads are moist and are making good contact.
During the massage, the skin feels a painful burning sensation, or the stimulation becomes weakened.	Adhesive gel pads are not adhering firmly to the skin or the gel pads are too dry.	Apply a few drops of water to adhesive surface of each pad and make sure the pads are pressed firmly to the skin during application.
Unit is on, but no sensation is felt in the pads.	Both pads are not firmly on the skin. Or there is a loose connection. Or the intensity level needs to be increased.	Make sure both pads are firmly pressed to the skin. Make sure all connections are secure from the unit to the wires and the wires to the pads. Increase the intensity level.
Adhesive gel pads do not stick to skin even after cleaning and moisturizing the gel pad.	Adhesive gel pads need to be replaced.	Replace the set of pads. Pads can be purchased from local distributor.
Unit does not turn on or the screen is dim.	Battery is low.	Recharge the battery.
Does not provide adequate pain relief.	Not being used long enough or individual is not lying down then massaging neck.	Use your massager 20-30 minutes at a time, 3-6 times per day. Lie on your back for neck support when you need relief for headaches or neck pain.

## Symbol Meaning



**1** Connect the wire to the re-usable adhesive pads.

**3** Place the pads on your skin.

**5** Turn on the device.

**2** Remove the plastic sheet over the pads.

**4** Plug the ports into the device.

**6** Select the Massage mode and Intensity.



UPPER BACK



DELTOID




ELBOW



HIP / THIGH



LOWER BACK



SCIATICA



SHOULDER / KNOTS




NECK



CARPAL TUNNEL



ANKLE / FOOT



KNEE



CALF